



HEALTHWORLD
INTERNATIONAL PHILS. CORP.

World Class QUALITY Need Not Be Expensive

Compact LED Light Source



Specifications:

- Voltage: AC95~240 50HZ~60HZ V
- Power: 50*1 W
- Life of Bulb: 500 h
- Illuminance: 3000000 lx
- Size: 100*190*280 mm
- Net Weight: 3.9 kg
- Package Size: 340*320*230 mm

Distributed By:

HEALTHWORLD INTERNATIONAL PHILIPPINES CORPORATION

📍 0421 Quirino Ave., La Huerta, Parañaque City, Metro Manila 1700 Philippines

☎ (+63) 2 8829-0175 | (+63) 919-065-7403 | (+63) 917-822-2938 (+63) 998-581-2446

✉ info@healthworldintl.com