



HEALTHWORLD
INTERNATIONAL PHILS. CORP.

World Class QUALITY Need Not Be Expensive

**PORTABLE HD
FLEXIBLE
RHINOLARYNGOSCOPE
WITH VIDEO PROCESSOR
(CHIP ON THE TIP)**



RP11A



FLEXIBLE VIDEO RHINOLARYNGOSCOPE



HD




HEALTHWORLD INTERNATIONAL PHILS. CORP.

World Class QUALITY Need Not Be Expensive

160g

Ultra-light handle, more stable grip, reduce the fatigue of the doctor, and make the examination smoother.



1280x720P

3.2mm high-definition video endoscope, ultra-thin insertion, full of details

A variety of connection methods are suitable for various diagnosis and treatment environments.



Computer



Pad



Processor

Optical system

■ Field of view	90°
■ Depth of field	3-50mm
■ Direction of view	0°

Insertion portion

■ Insertion tube OD	3.2mm
■ Angulation range	U 130° / D 130°
■ Working length	330mm±10%
■ Extension cable	USB、REMO
■ Max. resolution	1280x720P

Distributed By:

HEALTHWORLD INTERNATIONAL PHILIPPINES CORPORATION

📍 0421 Quirino Ave., La Huerta, Parañaque City, Metro Manila 1700 Philippines

☎ (+63) 2 8829-0175 | (+63) 919-065-7403 | (+63) 917-822-2938 (+63) 998-581-2446

✉ info@healthworldintl.com