




HEALTHWORLD
INTERNATIONAL PHILS. CORP.

World Class QUALITY Need Not Be Expensive

**PORTABLE HD
FLEXIBLE SCOPE
W/ VIDEO PROCESSOR
(CHIP ON THE TIP)**



HD



RF15A 
Portable High-Definition Digital Bronchoscope



HEALTHWORLD INTERNATIONAL PHILS. CORP.


World Class QUALITY Need Not Be Expensive

Versatility

It is suitable for trachea and difficult airway intubation, airway and double-lumen inspection. Alveolar lavage, biopsy sampling, pneumothorax closure, removal of tracheal foreign bodies and other bronchial treatments.

Applicability

Application departments: ICU, NICU, respiratory department, anesthesiology department, emergency department, etc., basically meet the clinical needs of each department



Optical system

■ Field of view	90°
■ Depth of field	3–50mm
■ Direction of view	0°

Insertion portion

■ Insertion tube OD	4.0mm/4.9mm/5.5mm
■ Angulation range	U 180° / D 130°
■ Inside diameter	1.2mm/2.2mm/2.6mm
■ Working length	600mm
■ Max. resolution	1280×720P

Distributed By:

HEALTHWORLD INTERNATIONAL PHILIPPINES CORPORATION

📍 0421 Quirino Ave., La Huerta, Parañaque City, Metro Manila 1700 Philippines

☎ (+63) 2 8829-0175 | (+63) 919-065-7403 | (+63) 917-822-2938 (+63) 998-581-2446

✉ info@healthworldintl.com