



# HEALTHWORLD INTERNATIONAL PHILS. CORP.

World Class *QUALITY* Need Not Be Expensive

## FDA PHILIPPINES-REGISTERED

**ENTITY**<sup>™</sup>  
Flexible Nasopharyngoscopes

## ENTity Flexible Fiberoptic Nasopharyngoscope

*Designed to empower ENT clinicians with the visibility & mobility they need*

The ENTity Series of Flexible Fiberoptic Nasopharyngoscopes can be used virtually anywhere you need to be. It is the only flexible endoscope that carries a guaranteed two (2) year warranty due to its reliability and durability. Our patented OptiLux integrated LED illumination system eliminates cumbersome external light boxes while providing exceptional visibility, making the endoscope completely portable. The OptiLux lighting system within the endoscope handle provides superior lighting targeted in the 5500K color temperature range for superior tissue color differentiation.



*Industry's only Nasopharyngoscope with Self-Contained Light Source & 2-Year Warranty*

### ENTITY SDXL

ENTity SDXL features a slimline 2.7mm insertion shaft making it the preferred choice for diagnostic evaluation of narrow upper airway passages in both adults and children.

|                      |                       |
|----------------------|-----------------------|
| Flexible Length      | 30cm (11.8")          |
| Overall Length       | 53cm (20.9")          |
| Shaft Diameter       | 2.7mm                 |
| Weight               | 11.6 oz               |
| End Tip Articulation | 110 degrees up / down |
| Depth of Field       | 5mm to 50mm           |
| Field of View        | 85 degrees            |
| Battery              | Lithium-Ion           |

### ENTITY XL

ENTity XL makes diagnostic examination of the nose and throat for viewing of the nasal passageways and larynx as easy to use as a stethoscope.

|                      |                       |
|----------------------|-----------------------|
| Flexible Length      | 30cm (11.8")          |
| Overall Length       | 53cm (20.9")          |
| Shaft Diameter       | 3.6mm                 |
| Weight               | 11.6 oz               |
| End Tip Articulation | 135 degrees up / down |
| Depth of Field       | 5mm to 50mm           |
| Field of View        | 70 degrees            |
| Battery              | Lithium-Ion           |

#### Better economics

- See more patients per day, increase revenue
- Reduced capital expenses and ongoing maintenance
- No lamp or light guide replacement costs

#### Portable, compact, all-in-one design

- No external light source required
- Fully sealed and completely soakable
- Ideal for office, bedside & hospital use

#### Distributed By:

**HEALTHWORLD INTERNATIONAL PHILIPPINES CORPORATION**

0421 Quirino Ave., La Huerta, Parañaque City, Metro Manila 1700 Philippines

(+63) 2 8829-0175 | (+63) 919-065-7403 | (+63) 917-822-2938 (+63) 998-581-2446

info@healthworldintl.com

Optim<sup>LLC</sup>



**OPTIM'S PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA**

Optim is ISO13485:2016 Certified